

European Data Format



Specs

Downloads

Companies



Home EDF > Specs > EDF+

Full specification of EDF+

The specification is also in the [original article](#) as published by [Elsevier](#): Bob Kemp and Jesus Olivan. European data format 'plus' (EDF+), an EDF alike standard format for the exchange of physiological data. Clinical Neurophysiology, 114 (2003): 1755-1761.

Acknowledgement

Many EDF users suggested to develop something like EDF+. We made a proposal in the summer and the specification was finalized in December 2002. We very much appreciate the constructive discussion with Stig Hanssen, Peter Jacobi, Kevin Menningen, Garðar Þorvarðsson, Thomas Penzel, Marco Roessen, Andreas Robinson and Alpo Värri, mainly on and around Yahoo's EDF users group.

Contents

1. [Introduction](#)
2. [The EDF+ protocol](#)
 - 2.1. [EDF+ compared to EDF](#)
 - 2.1.1. [Header](#)
 - 2.1.2. [Data records](#)
 - 2.1.3. [Additional specifications](#)
 - 2.2. [Annotations for text, time-keeping, events and stimuli](#)
 - 2.2.1. [The 'EDF Annotations' signal](#)
 - 2.2.2. [Time-stamped Annotations Lists \(TALs\)](#) in an 'EDF Annotations' signal
 - 2.2.3. [Annotations in a TAL](#)
 - 2.2.4. [Time keeping](#) of data records
 - 2.3. [Analysis results](#)
3. [Examples](#)
 - 3.1. [Auditory EP recording](#)
 - 3.2. [Sleep PSG and MSLT](#)
 - 3.3. [Sleep scoring](#)
 - 3.4. [Neurophysiology](#)
 - [Continuous EMG](#)
 - [F response](#)
 - [Motor Nerve Conduction](#)
 - [Somatosensory Evoked Potential](#)
 - [Visual Evoked Potential](#)
 - 3.5. [Intraoperative monitoring](#)
 - 3.6. [Montages in a routine EEG](#)
 - 3.7. [The Motor Nerve Conduction file](#)
 - [Header Record](#)
 - [Data Records](#)
 - [TALs](#)

1. Introduction

After its introduction in 1992, the [European Data Format \(EDF\)](#) became the standard for EEG and PSG (Sleep) recordings. During that time, several users pointed us at its limitations also for application in other fields such as myography, evoked potentials and cardiology. A major limitation was that EDF can only handle uninterrupted recordings. So we simply skipped that limitation but kept all other specifications of EDF intact. While maintaining EDF compatibility, we also standardized most of the labels and added a possibility to save annotations and analysis results. The result is EDF+ and can save most EEG, PSG, ECG, EMG, and Evoked Potential data that can not be saved into common hospital information systems.

Using EDF+, all signals, annotations and events that are recorded in one session using one recording system can be kept safely together in one file. EDF+ can also store events and annotations only, without any signals. This flexibility allows the user to choose an optimal mix. For instance, our sleep centre stores all on-line recorded data (signals, annotations) in one file, its hypnogram and apnea detections in another EDF+ file, the same sleep scorings but made by another technician in a third file. In Neurophysiology, the on-line obtained raw EEG traces with stimulus events from an EP investigation might be stored in one file and the averaged curves with detected latencies in a second file. In Cardiology, the raw ECG with annotations about the patients exercises can be in one file, the detected QRS parameters in another file.

EDF+ allows storage of several NON-CONTIGUOUS recordings into one file. This is the *only* incompatibility with EDF. All other features are EDF compatible. In fact, old EDF viewers still work and display EDF+ recordings as *if* they were continuous. Therefore, we recommend EDF+ files of *EEG* or *PSG studies* to be continuous if there are no good reasons for the opposite.

Because EDF+ is very close to EDF, and equally simple, EDF+ software can relatively easily be developed based on available EDF software.

2. The EDF+ protocol

Because EDF+ is based on EDF, you should first read the [EDF specs](#). Section 2.1 below describes how EDF+ differs from EDF. Section 2.2 describes how one of the EDF+ signals can be specially coded to store text annotations, time, events and stimuli.

EDF+ prescribes the following filename convention. Signals, recorded using the same technique and constant amplifier settings, can be stored in one file. Different techniques, or identical techniques but with different amplifier settings, must be stored in separate files. All EDF+ files must have .edf or .EDF as filename extension. See also section 2.3.

2.1. EDF+ compared to EDF

A standard EDF file consists of a header record followed by data records. The header record identifies the patient and specifies the technical characteristics of the recorded signals. The data records contain consecutive fixed-duration epochs of the recorded signals. A standard EDF+ file also consists of a header record followed by data records. The structure of these records is compatible to EDF but contains additional specifications. Note that for your EDF+ software to also be EDF compatible, it should support but *not rely* on these additional specifications.

2.1.1. The EDF+ header

The EDF+ header record identifies the patient and specifies the technical characteristics of the recorded signals according to the [EDF specs](#), except for the *first 'reserved' field (44 characters) which must start with 'EDF+C' if the recording is uninterrupted*, thus having contiguous data records, i.e. the starttime of each data record coincides with the end (starttime + duration) of the preceding one. In this case, the file is EDF compatible and the recording ends (number x duration) seconds after its startdate/time. The *'reserved' field must start with 'EDF+D' if the recording is interrupted*, so not all data records are contiguous. In both cases, the time must be kept in each data record as specified in [section 2.2.4](#).

The only incompatibility with EDF is, that signals may be recorded discontinuously. Therefore, we have decided that the EDF+ 'version' field must still read '0' 'like in EDF. In this way, old EDF viewers will still work and display EDF+ files (be they continuous or discontinuous) as continuous EDF files. EDF+ software will know the difference between continuous and discontinuous files from the mentioned 'reserved' field.

2.1.2. The EDF+ data records

A signal in an EDF data record is (a series of 2-byte samples, the subsequent samples representing subsequent integer values of that signal, sampled with equal time intervals. We will refer to this kind of signal as an '*ordinary signal*' from now on. EDF+ data records can and usually do) also contain ordinary signals. The EDF+ data records contain the ordinary signals according to the [EDF specs](#) (including the size limit of 61440), except for the fact that the data records may *unconditionally* be shorter than 1s, and *subsequent data records need not form a continuous recording*. However, as in EDF, data records that follow up in time must also follow up in the file. The samples of an ordinary signal must have equal sample intervals inside each data record, but the interval to the first sample of the next data record may be different.

For example, in a motor nerve conduction study with a number of stimuli, each data record would hold the ordinary signals corresponding to one stimulus. In this case, the duration of a data record corresponds to the "window size" in an ENMG/EP study.

Specifying a duration makes no sense if the EDF+ file does not contain any ordinary signals, for instance in a file that only contains manual sleep scores (like in both demo scoring files that come with [Polyman](#), and the [example in 3.3](#)). This is also true in the extreme case in which each ordinary signal only occupies one sample in each data record, while the file is discontinuous (EDF+D). *In those two cases, specify the 'duration of a data record' to be 0.*

2.1.3. Additional specifications in EDF+

1. In the header, use only printable US-ASCII characters with byte values 32..126.
2. The 'startdate' and 'starttime' fields in the header should contain only characters '0-9' and the period (.) as a separator, for example "02.08.51". In the 'startdate', use 1985 as a clipping date in order to avoid the Y2K problem. So, the years 1985-1999 must be represented by yy=85-99 and the years 2000-2084 by yy=00-84. After 2084, yy must be 'yy' and only item 4 of this paragraph defines the date.
3. The 'local patient identification' field must start with the subfields (subfields do not contain, but are separated by, spaces):
 - the code by which the patient is known in the hospital administration.
 - sex (English, so F or M).
 - birthdate in dd-MMM-yyyy format using the English 3-character abbreviations of the month in capitals. 02-AUG-1951 is OK, while 2-AUG-1951 is not.
 - the patients name.

Any space inside the hospital code or the name of the patient must be replaced by a different character, for instance an underscore. For instance, the 'local patient identification' field could start with: MCH-0234567 F 02-MAY-1951 Haagse_Harry. Subfields whose contents are unknown, not applicable or must be made anonymous are replaced by a single character 'X'. So, if everything is unknown then the 'local patient identification' field would start with: 'X X X X'.

- Additional subfields may follow the ones described here.
4. The 'local recording identification' field must start with the subfields (subfields do not contain, but are separated by, spaces):
 - The text 'Startdate'.
 - The starttime itself in dd-MMM-yyyy format using the English 3-character abbreviations of the month in capitals.
 - The hospital administration code of the investigation, i.e. EEG number or PSG number.
 - A code specifying the responsible investigator or technician.
 - A code specifying the used equipment.

Any space inside any of these codes must be replaced by a different character, for instance an underscore. The 'local recording identification' field could contain: Startdate 02-MAR-2002 PSG-1234/2002 NM Telemetry03. Subfields whose contents are unknown, not applicable or must be made anonymous are replaced by a single character 'X'. So, if everything is unknown then the 'local recording identification' field would start with: 'Startdate X X X X'. Additional subfields may follow the ones described here.

5. 'Digital maximum' must be larger than 'Digital minimum'. In case of a negative amplifier gain the corresponding 'Physical maximum' is smaller than the 'Physical minimum'. Check item 9 on how to apply the 'negativity upward' rule in Clinical Neurophysiology to the physical ordinary signal. 'Physical maximum' must differ from 'Physical minimum'. In case of uncalibrated signals, physical dimension is left empty (that is 8 spaces), while 'Physical maximum' and 'Physical minimum' still contain default values (this is to avoid division by 0 errors by some viewers).
6. Never use any digit grouping symbol in numbers. Never use a comma ',' for a or for a decimal separator. When a decimal separator is required, use a dot ('.').
7. The ordinary signal samples (2-byte two's complement integers) must be stored in 'little-endian' format, that is the least significant byte first. This is the default format in PC applications.
8. The 'starttime' should be local time at the patients location when the recording was started.
9. Use the standard texts and polarity rules at <http://www.edfplus.info/specs/edftexts.html>. These standard texts may in the future be extended with further texts, a.o. for Sleep scorings, ENG and various evoked potentials.

10. The 'number of data records' can only be -1 during recording. As soon as the file is closed, the correct number is high-pass and must be entered.

11. If filters (such as HighPass, LowPass or Notch) were applied to the ordinary signals then, preferably automatically, specify them like "HP:0.1Hz LP:75Hz N:50Hz" in the "prefiltering" field of the header. If the file contains an analysis result, the prefiltering field should mention the relevant analysis parameters.

12. The "transducertype" field should specify the applied sensor, such as "AgAgCl electrode" or "thermistor".

2.2. Annotations for text, time-keeping, events and stimuli

This section describes how one of the EDF+ signals can be specially coded to store text annotations, time, events and stimuli. In this way, annotations and events are kept in the same file as the signals that they refer to. The coding is EDF compatible in the sense that old EDF software would simply treat this 'EDF Annotations' signal as if it were a (strange-looking) ordinary signal.

2.2.1. The 'EDF Annotations' signal

EDF+ data records (and often do) contain ordinary signals. EDF+ introduces one other kind of signal, in which the values are annotations that can occur at any arbitrary point of time. This signal is identified by giving it (in the EDF+ header) the label 'EDF Annotations'. As in EDF, the '*nr of samples in each data record*' field in the header specifies how many 2-byte integers this 'EDF Annotations' signal occupies in each data record. But instead of storing 'ordinary signal' samples, those 2-byte integers are kept with characters. The character-bytes are stored byte-by-byte without changing their order. For instance, the text 'abc' is represented by successive bytes 97, 98 and 99 in the 'EDF Annotations' signal, though no annotations are to be kept, an EDF+ file must contain at least one 'EDF Annotations' signal in order to specify the starttime of each data record (see section 2.2.4). Of course, the label 'EDF Annotations' is not allowed for ordinary signals.

The 'EDF Annotations' signal only has meaningful header fields 'label' and 'nr of samples in each data record'. For the sake of EDF compatibility, the fields 'digital minimum' and 'digital maximum' must be filled with -32768 and 32767, respectively. The 'Physical maximum' and 'Physical minimum' fields must contain values that differ from each other. The other fields of this signal are filled with spaces.

2.2.2. Time-stamped Annotations Lists (TALs) in an 'EDF Annotations' signal

Text, time-keeping, events and stimuli are coded as text annotations in this 'EDF Annotations' signal. The annotations are listed in Time-stamped Annotations Lists (TALs) as follows. Each TAL starts with a time stamp OnsetDuration in which **00** and **00** are single bytes with value '21' and '20', respectively (unprintable ASCII characters) and Onset as well as Duration are coded using US-ASCII characters with byte values 43, 45, 46 and 48-57 (the '+', '-', '.', '0' and '9' characters, respectively). Onset must start with a '+' or a '-' character and specifies the amount of seconds by which the onset of the annotated event follows ('+') or precedes ('-') the startdate/time of the file, that is specified in the header. Duration must not contain any '+' or '-' and specifies the duration of the annotated event in seconds. If such a specification is not relevant, Duration can be skipped in which case its preceding **00** must also be skipped. Both Onset and Duration can contain a dot (.) but only if the fraction of a second is specified (up to arbitrary accuracy).

After the time stamp, a list of annotations all sharing the same Onset and Duration may follow. Each annotation is followed by a single **00** and may not contain any **00**. A **00**-byte (the unprintable ASCII character with byte value 0) follows after the last **00** of this TAL. So the TAL ends with a **00** followed by a **00**.

In each data record, the first TAL must start at the first byte of the 'EDF Annotations' signal'. Subsequent TALs in the same data record must follow immediately after the trailing **00** of the preceding TAL. A TAL, including its trailing **00**, may not overflow into another data record. Each event is annotated only once, even if its duration makes it extend into the time period of other data records. Unused bytes of the 'EDF Annotations' signal in the remainder of the data record are also filled with **00**-bytes. Additional 'EDF Annotations' signals may be defined according to the same specification.

For example, if the technician switches off the lights and closes the door 3 minutes after startdate/time, this can be stored as the 28-bytes TAL '+1800**00**lights off**00**+1800**00**Close door**00**' without the quotes. Alternatively, the two events can be stored as two separate shorter TALs '+180**00**lights off**00**+180**00**Close door**00**', also without the quotes. The TAL '+1800.2**00**25.5**00**Apnea**00**' codes a 25.5s apnea that begins 30 minutes and 0.2s after starttime.

2.2.3. Annotations in a TAL

The part between the next**00** is called one annotation. These annotations may only contain UCS characters (ISO 10646, the 'Universal Character Set', which is identical to the [Unicode](#) version 3+ character set) encoded by [UTF-8](#). This encoding is supported by the major operating systems, compilers and applications. The first 127 UCS characters are identical to those in US-ASCII and are encoded in the corresponding single byte values. US-ASCII characters that are represented by byte values 0-31 are allowed in the annotations only if explicitly prescribed by this EDF+ protocol. In order to enable multi-line texts and tables, US-ASCII characters that are represented by byte values 9 (TAB), 10 (LF) and 13 (CR) are allowed in the annotations. The first 65534 characters (the Basic Multilingual Plane: BMP) of the UCS contain virtually all characters used in any language in the world including Asian languages and UTF-8 encodes these in up to three byte-values. Remember that this encoding applies to the 'EDF Annotations' signal only: in the EDF+ header, only US-ASCII characters with byte values 32..126 are allowed.

In order to support automatic averaging and superimposition, identical events or stimuli that occur several times in one file must be coded each time by the same, unique annotation. Annotations (the part between **00** and the next **00**) of different events/stimuli (or the same stimulus at a different location) must differ from this unique annotation.

Annotations, for instance stimuli, that are related to information in only one particular data record, must be in that same data record. Even annotations describing events preceding the start of that data record, for instance a pre-interval stimulus must follow the time-keeping annotation.

2.2.4. Time keeping of data records

Because data records need not be contiguous, the starttime of each data record must be specified in another way. So, the *first annotation of the first* 'EDF Annotations' signal in each data record is empty, but its timestamp specifies how many seconds after the filestartdate/time that data record starts. So, if the first TAL in a data record reads '+567**00**', then that data record starts 567s after the startdate/time of the file. If the data records contain 'ordinary signals', then the analysis file should have startdate 03.08.99 and starttime 01.05.00. In this way it is clear that both files refer to one and the same time period in the patient's life. Some viewers (like PolyMan) are then capable of showing the two (or more) files time-synchronized on one screen. Because the analysis may reduce or increase the amount of data, the durations of analysis-file data records and recording-file data records may differ.

Apply suitable scaling factors in such a way that a large part of the available range of -32768 till 32767 for the values of the analysis results is used. If necessary, the scaling factor can be adapted to the dynamic range of the analysis result even after the analysis was done. Put these scaling factors in the header (digital and physical minimum and maximum) of the analysis file. If such scaling is really impossible because the useful dynamic range of the analysis result is too large, but only then, apply the [standardized logarithmic transformation](#) to store floating point values. Be aware that old EDF software is not aware of this transformation, and will show the analysis results on a logarithmic scale. So really try scaling first!

If a hypnogram is stored as an ordinary signal, sleep stages W,1,2,3,4,R,M should be coded in the data records as the integer numbers 0,1,2,3,4,5,6 respectively. Unscored epochs should be coded as the integer number 9. If a hypnogram is stored as annotations, use the [standardized](#) coding.

Automatically document the analysis principle and parameters in the Recording-id and, in case of ordinary signals, also in Label, Transducer type, Physical dimension and Prefiltering fields in the header of the analysis file.

3. Some examples

3.1. Auditory EP recording

The following is an example of annotations in the first two data records of an auditory EP recording. Each data record has two TALs, the first one includes the (mandatory) time-keeping annotation, the second one specifies a pre-interval stimulus.

	+0 00 Stimulus click 35dB both ears 00 Free text 00
	-0.065 00 Pre-stimulus beep 1000Hz 00
	+0.3 00 Stimulus click 35dB both ears 00
	+0.235 00 Pre-stimulus beep 1000Hz 00

In this example, averaging can be triggered by the unique texts "Stimulus click 35dB both ears" and/or "Pre-stimulus beep 1000Hz".

3.2. Sleep recording (PSG) with MSLT

A PSG that is followed by an MSLT can be stored in separate files. The PSG file, including lights-off and final wake-up annotations, is a continuous EDF+ file. The MSLT is a discontinuous EDF+ file which contains only the 20-minute periods in bed. Alternatively, the PSG and the MSLT can also be stored together into one single (discontinuous) file.

3.3. Sleep scoring

A 8-24hr sleep recording takes about 30-300MB when stored in EDF+. The recording can be analyzed manually, resulting in apnea's, leg movements and sleep stages. These results are kept in a separate EDF+ file (about 10-100kB) which, in this example, only contains one data record with one 'EDF Annotations' signal and no 'ordinary' signals. The table below shows the first half hour and the last few minutes in this data record. This patient fell asleep 9 minutes after switching off the light and had limb movements (Right and/or Left leg) and apneas after reaching sleep stage 2 and 3, respectively. If another technician also scores apnea's, leg movements or sleep stages, these scorings can be kept in another separate EDF+ file.

	+0 00 Recording starts 00
	+0 00 60 00 Sleep stage W 00
	+120 00 Lights off 00
	+660 00 300 00 Sleep stage N1 00
	+742 00 Turning from right side on back 00
	+960 00 180 00 Sleep stage N2 00
	+993.2 00 1.2 00 Limb movementR+L leg 00
	+1019.4 00 0.8 00 Limb movementR leg 00
	+1140 00 300 00 Sleep stage N3 00
	+1526.8 00 30.0 00 Obstructive apnea 00
	+1603.2 00 24.1 00 Obstructive apnea 00
	+1440 00 210 00 Sleep stage N2 00
	+1650 00 270 00 Sleep stage N3 00
	+1634 00 Turning from back on left side 00
	+1920 00 300 00 Sleep stage N2 00

	+3010 00 Lights on 00
	+30210 00 Recording ends 00

3.4. A large neurophysiological session

The example session includes the following investigations that are stored in separate EDF+ files. Note that the EDF+ processing software must implement any 'negativity upward' rule after reading the signal from the file. See further label and polarity rules at <http://www.edfplus.info/specs/edftexts.html>.

A continuous EMG is stored in a file with two signals: the raw EMG and the obligatory 'EDF Annotations' signal. This is a 'continuous' EDF+ file, so the recording might alternatively be stored as an EDF file. In case of a concentric needle electrode recording, a positivity at the centrally insulated wire relative to the cannula of the needle is stored as a positive value in the file.

An F response is also stored in a file with the raw EMG and the annotations. The duration of a data record equals the duration of the investigator's screen, i.e. of the 'window' (e.g. 50 ms). Each data record contains one single response. The annotations describe timing and any further characteristics of the stimulus. They can also describe measured distances and latencies.

A Motor Nerve Conduction Velocity, with one EMG channel is also stored in a file with the raw EMG and the annotations. The duration of a data record equals the duration of the 'window'. The curves from wrist and elbow stimulation are stored in the first and the second data record, respectively. The annotations describe timing and further characteristics of the stimulus. They can also describe measured distances and latencies. This file is described in full detail below.

A Somatosensory Evoked Potential with four recorded signals is stored in a file with five signals: the four raw SSEP signals and one 'EDF Annotations' signal. The duration of a data record equals the duration of the window (e.g. 100 ms). The annotations describe timing and characteristics of the stimulus. Another EDF+ file contains the 4-channel averaged responses (averages of odd and even sweeps are kept in different data records) and the 'EDF Annotations' signal which stores stimulus characteristics and measured latencies.

A Visual Evoked Potential is investigated by recording two sagittal EEG signals during checkerboard stimulation of the left and right field, respectively. The left- and right-stimulated averages are stored in two separate files. Left and right investigations are repeated once for checking reproducibility, so both files contain two data records. Each data record lasts 300ms and contains three signals: two EEG averages and the EDF Annotations. The annotations describe stimulus characteristics and measured latencies. The signal sampling starts 10ms before each stimulus, so the first two TALs in the 'left' file is 0.000**00** and 0.010**00**Stimulus checkerboard left**00**.

3.5. Intra-operative monitoring

Four (left and right) signals with alternating right and left stimulus are monitored. In this case, left and right stimulation affect correspondingly lateralized signals. So, the recording can be stored in two ('left' and 'right') EDF+ files. Each file contains 4 electrophysiological signals and 1 'EDF Annotations' signal. Alternatively, if sufficient amplifiers are available, the recording can be stored in one file containing 9 signals (4 'left' and 4 'right' electrophysiological signals and 1 'EDF Annotations' signal). In both cases, the response to each stimulus is stored in a separate data record and the 'EDF Annotations' signal specifies timing and characteristics of the stimulus (a.o. whether it is left or right).

3.6. Routine EEG

The 10/20 system electrodes (for instance F3, C3, T3, Cz and O1 and so on) are recorded against a common reference and saved as such in the EDF+ file. Therefore the montages can be made during review, such as F3-C3, T3-C3, C3-Cz and C3-O1. Because electrode locations are specified using [standard texts](#), re-montaging (i.e. re-referencing) EEG derivations can be done automatically. The Annotations signal contains events such as 'Eyes Closed' or 'Hyperventilation'.

3.7. The Motor Nerve Conduction file

A right Median Nerve conduction velocity is investigated by recording the right Abductor Pollicis Brevis while stimuli are given at wrist and elbow. The averaged signal and the corresponding annotations are stored in two data records: one for wrist and the other for elbow stimulation.

The header record contains	
8 ascii : version of this data format (0)	0
80 ascii : local patient identification	MCH-0234567 F 02-MAY-1951 Haagse_Harry
80 ascii : local recording identification	Startdate 02-MAR-2002 EMG5661 BK/JOP Sony. MNC R Median Nerve.
8 ascii : startdate of recording (dd.mm.yy)	17.04.01
8 ascii : starttime of recording (hh.mm.ss)	11:25.00
8 ascii : number of bytes in header record	768
44 ascii : reserved	EDF+D
8 ascii : number of data records (-1 if unknown)	2
8 ascii : duration of a data record, in seconds	0.050
4 ascii : number of signals (ns) in data record	2

ns * 16 ascii : ns * label		1st signal	2nd signal
ns * 80 ascii : ns * transducer type (e.g. AgAgCl electrode)		R APB	EDF Annotations
ns * 8 ascii : ns * physical dimension (e.g. uV)		mV	
ns * 8 ascii : ns * physical minimum (e.g. -500 or 34)		-100	-1
ns * 8 ascii : ns * physical maximum (e.g. 500 or 40)		100	1
ns * 8 ascii : ns * digital minimum (e.g. -2048)		-2048	-32768
ns * 8 ascii : ns * digital maximum (e.g. 2047)		2047	32767
ns * 80 ascii : ns * prefiltering (e.g. HP:0.1Hz LP:75Hz)		HP:3Hz LP:20kHz	
ns * 8 ascii : ns * nr of samples in each data record		1000	60
ns * 32 ascii : ns * reserved			

Each data record contains:
1000 * 2-byte integer : R APB samples
60 * 2-byte integer : EDF Annotations

The EDF Annotations signal in the 1st data record contains one TAL and is then filled out with **00**-bytes until the end. The TAL is:

+0 00 Stimulus right wrist 0.2ms x 8.2mA at 6.5cm from recording site 00 Response 7.2mV at 3.8ms 00
--

The EDF Annotations signal in the 2nd data record also contains one TAL:

+10 00 Stimulus right elbow 0.2ms x 15.3mA at 28.5cm from recording site 00 Response 7.2mV at 7.8ms (55.0ms) 00
--

In this example, the TALs take less than 100 characters per data record, so the header reserves 120 characters (60 'samples') for the EDF Annotations signal.

If desired, an internal structure of the annotations inside a TAL can be chosen, this is not obligatory, usually increases size, does not improve exchange to other systems, but can be useful locally. For example, the last TAL can be coded as four separate TALs and the annotation inside each TAL can be coded in XML, as follows:

+10 00 Stimulus_elbow 00
+10 00 EDF_XMLnote 00
<Stimulus_elbow><duration unit="ms">0.2</duration> <intensity mode="current" unit="mA">15.3</intensity> <position-right elbow>position> <distance mode="stimulus to recording" unit="cm">28.5</distance> </Stimulus_elbow> <EDF_XMLnote> 00
+10 00 <EDF_XMLnote> <measurements> <latency unit="ms">7.8</latency> <amplitude mode="baseline to peak" unit="mV">7.2</amplitude> <velocity mode="segmental" unit="m/s">55.0</velocity> </measurements> <EDF_XMLnote> 00